

SleepViewSM Direct

Home sleep testing fulfillment service designed for sleep centers.



Reduces operational burden

Enables providers to offer home sleep testing without managing logistics, inventory, or patient support.

Provides end-to-end service support

Patient scheduling and study support. Data quality review and notifications to lab staff.

Scalable testing capacity

Enable growth in testing programs and increase in study volume. Free up staff time to focus on scoring and patient care.

4,500+ Active users

97% Study success rate

95% patient satisfaction

A fully managed testing workflow:



1. Order

Provider determines patient eligibility and submits order



2. Ship

SleepView Direct schedules and ships the device to patients



3. Test

Patient completes the test at home and returns the device



4. Review

Data is processed and made available for scoring and interpretation

Expand your home sleep testing capabilities

Visit our website: beacon.bio or reach out to our team: dx.sales@beacon.bio

SleepView[®] Monitor

Home sleep testing system designed to support clinical evaluation with reliable data collection and analysis. SleepView[®] is a small and lightweight home sleep testing monitor that meets AASM recommended guidelines and is designed for patient self-administration.



Device specifications

Device type:	Type III (also meets Type IV requirements)
Channel count:	8
Dimensions:	3" x 2.6" x 0.7"
Weight:	2 oz.
Power:	1 AAA battery (powers 1 sleep study)
Memory:	1 GB internal memory (supports multiple-night recording)

Physiological signals collected

- Heart rate
- Pulse oximetry
- Respiratory airflow (cannula, CPAP compatible)
- Respiratory effort (RIP belt)
- Snore
- Body position
- Integrated actigraphy (total sleep time)

SleepView Web Portal

A complete home sleep testing workflow

- Centralized management of patient sleep studies and reports in a HIPAA-compliant database
- Online review and manual scoring of raw data
- Customizable report generation with digital signature
- Actigraphy-based sleep/wake analysis
- Automated respiratory event scoring
- Software upgrades and support included

